

Product data sheet

www.dinolite-uk.com

Dino-Lite MS12C stand

Variable view range caddy with two variable view range adapters, transparent. The MS12C is not compatible with polarizer models or models in the AD Pro 2 range.

The MS12C includes a base unit with two insert adapters providing three "hands free" fixed working distance and magnification options. Ideal for collectors, hobbyists, graphic artists as well as at home or in the classroom. The Dino-Lite digital microscopes simply slide into the base or one of the adapters. The clear acrylic base design allows easy accessibility of tools and light to the observed object.



Technical data

Model	MS12C stand
Unit dimensions (WxLxH)	8.26 x 8.26 x 7.62cm
Other	does not work with Polarizer and/or AD Pro 2 range models