

Product data sheet

www.dinolite-uk.com

Dino-Lite AM4113T5

The AM4113T5 has a fixed 500x magnification for high power applications. It is the ideal tool for examining very small objects under high magnification.



Technical data

Model	AM4113T5
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	500x fixed focus magnification
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	√ (by software)
LED color	White
Material	Composite
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√ (line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	-
Compatibility	Windows XP, Vista, 7 and recent Mac OS
Software (included)	Windows: DinoCapture, Mac: DinoXcope
Other	For examining very small objects when high magnification is needed. Scientific, biological or specialized industrial use.