

## Product data sheet

[www.dinolite-uk.com](http://www.dinolite-uk.com)

### Dino-Lite AM2011

The Dino-Lite AM-2011 is the basic model USB microscope with magnification of 10x -60x and 200x. The AM2011 features improved color rendering, compared to its predecessor. It is equipped with a high-speed USB 2.0 interface.



### Technical data

Model	AM2011
Resolution	640x480pixels (VGA)
Magnification	10x – 60x & 200x
Interface	USB 2.0
Nr. of LED's	4
LED on/off switchable	-
LED color	White
Material	Composite
Microtouch	-
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second
Measurement	-
Calibration	-
Filter/Diffuser	-
Compatibility	Windows XP, Vista, 7 and Mac OS 10.4 and newer
Software (Included)	Windows: DinoCapture 2.0, Mac: DinoXCope
Special	Affordable solution for home or school use. Suitable for hobbies & collections (coins, stamps, watches, jewels, etc.)