

Dino-Lite MS23B desktop clamp

The MS-23B articulating desktop clamp allows you to clamp a stand onto a tabletop or attach to any surface up to 3cm thick. The arm has a reach of 35cm allowing good distance and freedom to move the arm into flexible positions.



Technical data

Model	MS23B
Gooseneck length	35 cm long reach
Gooseneck thickness	1.6 cm diameter
Requirements	Can be clamped on surfaces of up to 3 cm thick
Unit weight	275g
Box dimensions	19x19x6 cm
Package weight	350g