

Product data sheet

www.dinolite-uk.com

Dino-Lite MS52BA2 stand

Heavy duty table base for the MS52B.

The MS52BA2 is a combination of the rigid MS52B articulating flex-arm and a desktop base rack that holds the MS52B on top of a table. The desktop base rack is designed for situations where the standard clamp cannot be used to secure the MS52B on a stable surface.

The flex-arm is moved and adjusted with tightening knobs and clamps for positioning the arm to the area you wish to examine. Two of the flex-arm axis can be rotated 360 degrees and extended for further reach.

The MS52BA2 flex-arm is manufactured from durable, lightweight recycled aluminum and has the ability to be folded flat, which makes the MS52BA2 stand portable and ready for travel.



Technical data

Model	MS52BA2
Unit dimensions (WxLxH)	36 x 15,3 x 5 cm
Unit weight	appr. 2270 g
Horizontal & vertical work area	621 - 921 mm
Top swivel	360 degrees
Other	-