

## Product data sheet

[www.dinolite-uk.com](http://www.dinolite-uk.com)

### Dino-Lite AM4113TL-M40

The AM4113TL-M40 is a special version of the long working distance series with a lower magnification rate of between 5 and 40 times.



### Technical data

|                       |   |
|-----------------------|---|
| <b>Model</b>          | <b>AM4113TL-M40</b>                       |
| Resolution            | 1,3 megapixel (1280x1024 pixels)          |
| Magnification         | 5x - 40x                                  |
| Interface             | USB 2.0, High speed USB transmission      |
| Nr. of LED's          | 8   |
| LED on/off switchable | √ (by software)                           |
| LED color             | White                                     |
| Material              | Composite                                 |
| Microtouch            | √   |
| Output                | Picture, video and time-lapsed video      |
| Video Frame Rate      | up to 30 frames-per-second, 15fps at 1.3M |
| Measurement           | √   |
| Calibration           | √   |
| Filter/Diffuser       | -   |
| Compatibility         | Windows XP, Vista, 7 and recent Mac OS    |
| Software (included)   | Windows: DinoCapture, Mac: DinoXcope      |