

Product data sheet

www.dinolite-uk.com

Dino-Lite AD4013MTL

The AD series have a detachable cap and is delivered with both a closed and an open exchangeable cap. This MTL model has a long working distance and aluminum alloy housing. The enhanced color sensor provides superior color rendering and improved low light imaging.



Technical data

Model	AD4013MTL
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	20x – 90x
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	√ (by software)
LED color	White
Material	Full aluminum alloy body
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√(line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	-
Compatibility	Windows XP, Vista, 7 and recent Mac OS
Software (Included)	Windows: DinoCapture, Mac: DinoXcope
Other	delivered with both a closed and an open exchangeable cap