

Product data sheet

www.dinolite-uk.com

Dino-Lite AM3013T

The AM3013T has both the measurement functions and the MicroTouch button to capture pictures directly from the Dino-Lite. The improved image sensor provides a high quality image with very natural colors and it also offers a high refresh rate for better real time images. It is an affordable solution for industrial, scientific and business use.



Technical data

Model	AM3013T
Resolution	640x480pixels (VGA)
Magnification	10x – 60x & 200x
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	√ (by software)
LED color	White
Material	Composite
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second
Measurement	√(line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	-
Compatibility	Windows XP, Vista, 7 and recent Mac OS
Software (Included)	Windows: DinoCapture, Mac: DinoXcope
Special	Affordable solution for industrial, scientific and business use