

Product data sheet

www.dinolite-uk.com

Dino-Lite AM4013MTL

Dino-Lite X Pro with aluminum alloy housing and long working distance and Microtouch. The aluminum housing offers enhanced protection and durability and has the best look and feel. The enhanced color sensor provides superior color rendering and improved low light imaging.



Technical data

Model	AM4013MTL
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	20x – 90x
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	√ (by software)
LED color	White
Material	Full aluminum alloy body
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√ (line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	-
Compatibility	Windows XP, Vista, 7 and recent Mac OS
Software (included)	Windows: DinoCapture, Mac: DinoXcope
Special	best model for repair, rework or assembly tasks